



Download on the
App Store



GET IT ON
Google Play

Scan to download Signia App

Scan to Contact



Scan to Email



The Information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice. The required features should therefore be specified in each individual case at the time of conclusion of the respective contract.

Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc.

“The Bluetooth” word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by WSAUD A/S is under license. Other trademarks and trade names are those of their respective owners. The “Qi” symbol is a trademark of the Wireless Power Consortium.”

Sivantos India Pvt. Ltd.

Sattva Knowledge Court,
Ground Floor, Wing ‘A’ - Office 1,
Plot No. 9, Doddenakkundi I Phase Industrial Area,
Doddenakkundi Village, K.R. Puram, Hoodi,
Bengaluru, Karnataka – 560048

Order ID: D00371498

Find us on:



www.signia.net/en-in



www.youtube.com/signiahearing/



www.facebook.com/signiaindia/



www.instagram.com/signiahearingindia/

Retail Price List

Effective from 17th March 2026



Be
Brilliant™

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, but performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath.

With Signia, you can Be Brilliant.

CONTENTS

1) Integrated Xperience (IX)	4
2) Augmented Xperience (AX)	18
3) Orion C&G	24
4) Sirion Connect	26
5) primax (px)	28
6) Other Hearing Aids	29
7) Accessories	30

Terms & Conditions

1. Signia Hearing Aids carry warranty as indicated in the price list.
2. Kit includes 2 units of hearing aids and a charger.
3. Accessories (Chargers, miniReceiver, Remote Control, Streamers, CROS / BiCROS) carry 1 year warranty.
4. Power adaptor for some accessories will have to be procured separately.
5. Custom hearing aids prices include shelling charges.
6. HP miniReceiver customization charges are included in the RIC pricing.
7. RIC / BTE instruments will be supplied in standard beige colour. Other colour housing options for RIC / BTE models are subject to a delivery lead time of 3 weeks.
8. Prices of color conversion kit may vary depending on hearing aid model.
9. Hearing Aids are exempted from GST.
10. GST applicable on accessories. Retail prices of accessories are inclusive of GST.
11. The pricing and product specifications are subject to change without prior notice.
12. Chargers for selective models will have to be purchased separately.

Sivantos India recommends all its valued consumers to seek a hearing care professional's opinion on the best suited hearing instrument. Online purchase may not necessarily give you the desired results.

Feature overview

Pure Charge&Go IX family, Styletto IX, Active Pro IX,
Motion Charge&Go IX family, Insio ITC[†]/ITE IX

Features	Description	Performance levels						
		Styletto IX, Pure Charge&Go IX family, Motion Charge&Go IX family, Insio ITC [†] /ITE IX					Active Pro IX	
		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
RealTime Conversation Enhancement	Provides uninterrupted engagement in noisy group conversations by seamlessly tracking the locations of all conversation partners and enhancing their speech. ¹	●●●●●	●●●●○	●●○○○			●●●●●	
Conversation Dynamics Analyzer	Analyzes conversation dynamics including turn-taking and timing, locations of pure noise without speech, and levels of background noise.	✓					✓	
Conversation Source Analyzer	Constantly scouts the acoustic environment to identify number of talkers and their location.	✓	✓	✓			✓	
Dynamic Focus Streams	Isolate and track speech independently.	3	1				3	
Conversation Booster	Dynamically enhances individual speakers.	●●●●●	●●●●○	●●○○○			●●●●●	
Augmented Focus	Separates the focus sound (speech) from surrounding sounds, creating a strong contrast. ¹	✓	✓	✓	✓	✓	✓	✓
Focus Stream	Enhances speech signals from the front hemisphere.	1	1	1	1	1	1	1
Surrounding Stream	Processes environmental sounds to preserve environmental awareness without masking speech.	1	1	1	1	1	1	1
Dynamic Soundscape Processing 3.0	Integrates inputs for precise hearing aid optimization, including (when available): acoustic sensors, motion sensors, OVP 2.0, noise detection, SNR analysis, and RealTime Conversation Enhancement.	●●●●●	●●●●○	●●●○○	●●○○○	●●○○○	●●●●●	●●○○○
Spatial SpeechFocus	Steers directional microphone beam to the dominant speech source (left, right, or behind the wearer automatically in a car situation or in a dedicated "Strolling" program). ¹	✓					✓	
SpeechFocus	Depending on the direction of the dominant speech source, Speech Focus steers the directional beam behind the wearer, or engages a surround-sound mode for speech from either side.	✓	✓				✓	
Own Voice Processing 2.0 (OVP™ 2.0) ²	Uses real-time recognition of the wearer's voice to deliver a natural own-voice impression. ¹	✓	✓	✓				
Motion sensor ³	Detects whether the wearer is walking, running, or sitting still to guide how sound is optimized.	✓	✓	✓			✓	
Binaural Directionality	Provides scientifically proven enhanced speech understanding in noise. ¹	✓	✓	✓			✓	
Adaptive Directionality	Improves the signal to noise ratio in noisy environments by reducing dominant noise sources from the back and sides. The polar plot is adjusted continuously in multiple frequency channels.	✓	✓	✓	✓		✓	
Automatic Directionality	Improves the signal to noise ratio in noisy environments by automatically and seamlessly applying directionality in individual processing channels.	✓	✓	✓	✓	✓	✓	✓
TruEar 360	Reproduces benefits of the natural pinna effect by highlighting the difference between front and back sources of sound.	✓	✓	✓	✓	✓	✓	✓

1)

Available for bilateral fittings with e2e wireless link

2)

Not available in Insio ITC/ITE IX and Active Pro IX

3)

Not available in Insio ITC/ITE IX

†)

Two-microphone models only

Features	Description	Performance levels						
		Styletto IX, Pure Charge&Go IX family, Motion Charge&Go IX family, Insio ITC [†] /ITE IX					Active Pro IX	
		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
SpeechStabilizer	Preserves the natural envelope of speech in any situation, including very noisy situations, which preserves acoustic cues that are important for speech understanding.	✓	✓	✓	✓	✓	✓	✓
Speech and Noise Management	Reduces listening effort using sophisticated sound processing.	✓	✓	✓	✓	✓	✓	✓
e2e wireless 4.0	Controls coupling, synchronization, CROS/BiCROS, and a range of sound features. ¹	✓	✓	✓	✓	✓	✓	✓
Listening Comfort								
Extended bandwidth ⁴	12 kHz available bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓					✓	
Auto EchoShield	Automatically reduces the effect of reverberations and sudden loud noises in reverberant environments.	✓					✓	
Reverberant Room Program	Dedicated program for reverberant environments.	✓	✓				✓	
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	●●●●●	●●●●●	●●●●●	●●●○○	●○○○○	●●●●●	●○○○○
eWindScreen™	Reduces the interference of wind noise when outdoors.	●●●●●	●●●●●	●●●○○	●○○○○	●○○○○	●●●●●	●○○○○
Extended dynamic range	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓	✓	✓
Feedback cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓	✓	✓
Wearer Empowerment								
Spatial Configurator	Allows the wearer to adjust the direction of the microphone beam via the Signia App. ¹	✓	✓				✓	
My WellBeing ²	Accessed via the Signia App, wearers can get insights about how much they are using their hearing aids, their level of social interaction, their daily step count, and their overall activity level.	✓	✓	✓			✓	
Signia Assistant ²	Offers wearers instant 24/7 troubleshooting and customized sound adjustments via an intuitive AI-enabled chat. Provides detailed information in Connexx to address each client's fine-tuning needs more precisely.	✓	✓	✓	✓	✓	✓	✓
Signia App (iOS and Android) ³	Easy-to-use smartphone app enables remote control adjustments, audio streaming, access to My WellBeing, and support via Signia Assistant and TeleCare.	✓	✓	✓	✓	✓	✓	✓
Adaptive Streaming Volume ²	Speech intelligibility weighted adaptation of the streaming volume depending on the level of ambient noise.	✓	✓	✓	✓	✓	✓	✓
TeleCare ³	Provides remote services including the TeleCare Dashboard, Remote Fine Tuning, and CareChat (convenient text messages, phone calls, and video calls).	✓	✓	✓	✓	✓	✓	✓

1)

Available for bilateral fittings with e2e wireless link

2)

Not available in Insio ITC/ITE IX

3)

Reduced features available for Insio ITC/ITE IX

4)

Not available in Motion Charge&Go SP IX

†)

Two-microphone models only

Feature overview

Pure Charge&Go IX family, Styletto IX, Active Pro IX,
Motion Charge&Go IX family, Insio ITC[†]/ITE IX

Features	Description	Performance levels						
		Styletto IX, Pure Charge&Go IX family, Motion Charge&Go IX family, Insio ITC [†] /ITE IX					Active Pro IX	
		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
Connectivity								
Direct audio streaming ²	Bluetooth Classic connectivity available in Pure Charge&Go BCT IX. MFi, LE Audio, and ASHA connectivity available in Pure Charge&Go IX / T IX, Styletto IX, Motion Charge&Go IX BTE family and Active Pro IX.	✓	✓	✓	✓	✓	✓	✓
HandsFree ²	Enables hands-free calling. The wearer's voice is picked up by the hearing aid microphone and delivered in high quality to the conversation partner, while the call audio is streamed directly to the wearer's hearing aids. Available for connections via Bluetooth Classic, MFi, or LE Audio.	✓	✓	✓	✓	✓	✓	✓
CallControl ²	Calls can be answered and ended by a press of the hearing aid button or by tapping the hearing aid, depending on hearing aid model. Available for connections via Bluetooth Classic, MFi, or LE Audio.	✓	✓	✓	✓	✓	✓	✓
Telecoil	Available in Pure Charge&Go BCT IX, Pure Charge&Go T IX, the Motion Charge&Go IX BTE family, and as an option in Insio ITC/ITE IX. Enables direct connection to audio loop systems and compatible telephone signals.	✓	✓	✓	✓	✓		
StreamLine Mic compatibility ¹	Stream high quality stereo sound directly from any Bluetooth device (PCs, older Apple devices, Android devices without ASHA or LE Audio, etc.). Can also be used as a remote mic.	✓	✓	✓	✓	✓	✓	✓
TV Sound [®] compatibility ²	Stream high quality stereo sound directly from TV to hearing aids via Bluetooth.	✓	✓	✓	✓	✓	✓	✓
miniPocket compatibility	Optional keychain remote control.	✓	✓	✓	✓	✓	✓	✓
Tinnitus								
Ocean Waves Therapy Signal (presets)	Therapy signals that mimic the sound of the sea to provide a soothing listening experience and reduce the impact of tinnitus.	4	4	4	4		4	
Static Therapy Signal (channels / presets)	Soft, gentle therapy signals that reduce the impact of tinnitus.	20/5	18/5	16/5	12/5		20/5	
Notch Therapy ³	Signia's unique, patented tinnitus therapy which is proven especially effective against tonal tinnitus.	✓	✓	✓	✓		✓	

1) Not compatible with Insio ITC/ITE IX or Pure Charge&Go BCT IX

2) Not available in Insio ITC/ITE IX

3) Not available in Motion Charge&Go SP IX

†) Two-microphone models only

Features	Description	Performance levels						
		Styletto IX, Pure Charge&Go IX family, Motion Charge&Go IX family, Insio ITC [†] /ITE IX					Active Pro IX	
		7 IX	5 IX	3 IX	2 IX	1 IX	Active Pro	Active
Fitting Optimization								
Processing Channels / Gain & MPO handles	Number of processing channels/handles in Connexx	48/20 ¹	36/18	32/16	24/12	16/8	48/20	16/8
HD Music	Preset optimized music programs for enhanced enjoyment of non-streamed music.	3	3	1	1		3	
Hearing programs	Number of programs configurable in Connexx.	6	6	6	4	4	6	4
CROS/BiCROS compatibility	A CROS device worn on an ear with unaidable hearing loss transmits sound to a hearing aid worn on the better ear.	✓	✓	✓	✓	✓	✓	✓
Frequency Compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓	✓	✓
Smart Optimizer and Data Logging	Tailored fine-tuning recommendations based on logging and analysis of the wearer's hearing aid use.	✓	✓	✓	✓	✓	✓	✓
Acclimatization Manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓	✓	✓
InSituGram	Allows quick and convenient measurement of hearing thresholds with the hearing aids in place.	✓	✓	✓	✓	✓	✓	✓
AutoFit [™]	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets, enabling significantly faster real-ear verification.	✓	✓	✓	✓	✓	✓	✓

1) Motion Charge&Go SP IX has 36/18 channels/handles

†) Two-microphone models only

Feature overview

Silk Charge&Go IX, Insio Charge&Go CIC IX, Insio IIC/CIC/ITC[†] IX

Features	Description	Performance levels				
Conversation management		7 IX	5 IX	3 IX	2 IX	1 IX
Dynamic Soundscape Processing 3.0	Integrates inputs for precise hearing aid optimization, including (when available): acoustic sensors, noise detection and SNR analysis.	●●●●●	●●●●○	●●●○○	●●○○○	●●○○○
Binaural OneMic Directionality 2.0	Microphones from each hearing aid work together to create a directional listening pattern preferentially amplifying speech from the front to support hearing in noisy environments. ¹⁾	✓	✓	✓		
SpeechStabilizer	Preserves the natural envelope of speech in any situation, including very noisy situations, which preserves acoustic cues that are important for speech understanding.	✓	✓	✓	✓	✓
Speech and Noise Management	Reduces listening effort using sophisticated sound processing.	✓	✓	✓	✓	✓
e2e wireless 4.0	Controls coupling, synchronization, CROS/BiCROS, and a range of sound features. ¹⁾	✓	✓	✓	✓	✓
Listening Comfort						
Extended bandwidth	12 kHz available bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓				
Auto EchoShield	Automatically reduces the effect of reverberations and sudden loud noises in reverberant environments.	✓				
Reverberant Room Program	Dedicated program for reverberant environments.	✓	✓			
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	●●●●●	●●●●●	●●●●●	●●●○○	●○○○○
eWindScreen™	Reduces the interference of wind noise when outdoors.	●●●●●	●●●●●	●●●○○	●○○○○	●○○○○
Extended dynamic range	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓
Feedback cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓
Wearer Empowerment & Connectivity						
Signia App (iOS and Android)	Easy-to-use smartphone app enables remote control adjustments and TeleCare.	✓	✓	✓	✓	✓
TeleCare	Provides remote services including the TeleCare Dashboard and CareChat (convenient text messages, phone calls, and video calls).	✓	✓	✓	✓	✓
miniPocket compatibility	Optional keychain remote control.	✓	✓	✓	✓	✓
Telecoil	Available as an option for Insio ITC IX. Enables direct connection to audio loop systems and compatible telephone signals.	✓	✓	✓	✓	✓

Features	Description	Performance levels				
Tinnitus		7 IX	5 IX	3 IX	2 IX	1 IX
Ocean Waves Therapy Signal (presets)	Therapy signals that mimic the sound of the sea to provide a soothing listening experience and reduce the impact of tinnitus.	4	4	4	4	
Static Therapy Signal (channels / presets)	Soft, gentle therapy signals that reduce the impact of tinnitus.	20/5	18/5	16/5	12/5	
Notch Therapy	Signia's unique, patented tinnitus therapy which is proven especially effective against tonal tinnitus.	✓	✓	✓	✓	
Fitting Optimization						
Processing Channels / Gain & MPO handles	Number of processing channels/handles in Connexx	48/20	36/18	32/16	24/12	16/8
HD Music	Preset optimized music programs for enhanced enjoyment of non-streamed music.	3	3	1	1	
Hearing programs	Number of programs configurable in Connexx.	6	6	6	4	4
CROS/BiCROS compatibility	A CROS device worn on an ear with unaidable hearing loss transmits sound to a hearing aid worn on the better ear. ²⁾	✓	✓	✓	✓	✓
Frequency Compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓
Smart Optimizer and Data Logging	Tailored fine-tuning recommendations based on logging and analysis of the wearer's hearing aid use.	✓	✓	✓	✓	✓
Acclimatization Manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓
InSituGram	Allows quick and convenient measurement of hearing thresholds with the hearing aids in place.	✓	✓	✓	✓	✓
AutoFit™	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets, enabling significantly faster real-ear verification.	✓	✓	✓	✓	✓

1) Available for bilateral fittings with e2e wireless link

†) One-microphone models

1) Available for bilateral fittings with e2e wireless link

2) Requires a hearing aid with e2e wireless link

RealTime Conversation Enhancement

The science of keeping up with conversations

Signia Integrated Xperience's revolutionary RealTime Conversation Enhancement analyzes, augments, and adapts to the dynamic flow of conversations. This lets wearers enjoy a seamless and enriched hearing experience.



Analyze

Processes 192,000 data points every second to identify speakers and analyze conversation dynamics, authentically reflecting real-life conversations.




Augment

Uses data gathered to control multiple focus streams that lock on to and dynamically highlight each speaker's voice amidst the noise. Meanwhile, Augmented Focus technology[2] processes background sound to ensure a comfortable, immersive conversational experience without distractions.




Adapt

Our patented multi-stream architecture updates 1000 times per second, merging multiple focus streams into a dynamic auditory space that adapts to evolving conversations. This allows for real-time tracking and enhancement of multiple moving speakers, ensuring wearers can fully engage in group discussions.



IX Power the conversation
with Signia Integrated Xperience



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
INTEGRATED XPERIENCE					

KIT STYLETTO IX

Kit Styletto 7IX	Slim-RIC	48	10 - 110	4	7,24,990
Kit Styletto 5IX	Slim-RIC	36	10 - 110	4	4,99,990
Kit Styletto 3IX	Slim-RIC	32	10 - 110	2	3,49,990
Kit Styletto 2IX	Slim-RIC	24	10 - 110	2	2,09,990
Kit Styletto 1IX	Slim-RIC	16	10 - 110	2	1,49,990

KIT PURE C&G BCT IX

Kit Pure C&G BCT 7IX	RIC	48	10 - 115	4	7,99,990
Kit Pure C&G BCT 5IX	RIC	36	10 - 115	4	5,29,990
Kit Pure C&G BCT 3IX	RIC	32	10 - 115	2	3,69,990
Kit Pure C&G BCT 2IX	RIC	24	10 - 115	2	2,09,990
Kit Pure C&G BCT 1IX	RIC	16	10 - 115	2	1,49,990

PURE C&G BCT IX

Pure C&G BCT 7IX	RIC	48	10 - 115	4	4,09,990
Pure C&G BCT 5IX	RIC	36	10 - 115	4	2,79,990
Pure C&G BCT 3IX	RIC	32	10 - 115	2	1,89,990
Pure C&G BCT 2IX	RIC	24	10 - 115	2	1,09,990
Pure C&G BCT 1IX	RIC	16	10 - 115	2	84,990

KIT PURE C&G IX / KIT PURE C&G T IX

Kit Pure C&G 7IX / T 7IX	RIC	48	10 - 115	4	7,69,990
Kit Pure C&G 5IX / T 5IX	RIC	36	10 - 115	4	4,79,990
Kit Pure C&G 3IX / T 3IX	RIC	32	10 - 115	2	3,46,990
Kit Pure C&G 2IX / T 2IX	RIC	24	10 - 115	2	2,04,990
Kit Pure C&G 1IX / T 1IX	RIC	16	10 - 115	2	1,44,990

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.

IX Power the conversation
with Signia Integrated Xperience



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
INTEGRATED XPERIENCE					

PURE C&G IX / PURE C&G T IX

Pure C&G 7IX / T 7IX	RIC	48	10 - 115	4	3,89,990
Pure C&G 5IX / T 5IX	RIC	36	10 - 115	4	2,49,990
Pure C&G 3IX / T 3IX	RIC	32	10 - 115	2	1,79,990
Pure C&G 2IX / T 2IX	RIC	24	10 - 115	2	1,07,990
Pure C&G 1IX / T 1IX	RIC	16	10 - 115	2	77,990

KIT MOTION C&G M/P/SP IX

Kit Motion C&G M/P/SP 7IX	BTE	48/48/36	10 - 120	4	7,09,990
Kit Motion C&G M/P/SP 5IX	BTE	36/36/36	10 - 120	4	4,45,990
Kit Motion C&G M/P/SP 3IX	BTE	32/32/32	10 - 120	2	2,39,990
Kit Motion C&G M/P/SP 2IX	BTE	24/24/24	10 - 120	2	1,79,990
Kit Motion C&G M/P/SP 1IX	BTE	16/16/16	10 - 120	2	1,34,990

MOTION C&G M/P/SP IX

Motion C&G M/P/SP 7IX	BTE	48/48/36	10 - 120	4	3,64,990
Motion C&G M/P/SP 5IX	BTE	36/36/36	10 - 120	4	2,32,990
Motion C&G M/P/SP 3IX	BTE	32/32/32	10 - 120	2	1,29,990
Motion C&G M/P/SP 2IX	BTE	24/24/24	10 - 120	2	94,990
Motion C&G M/P/SP 1IX	BTE	16/16/16	10 - 120	2	72,990

KIT ACTIVE IX

Kit Active Pro IX	Instant Fit ITE	48	10 - 80	4	3,24,990
Kit Active IX	Instant Fit ITE	16	10 - 80	2	1,59,990

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.

IX Power the conversation
with Signia Integrated Xperience



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
INTEGRATED XPERIENCE					
KIT SILK C&G IX					
Kit Silk C&G 7IX	Instant Fit CIC	48	10 - 80	4	7,24,990
Kit Silk C&G 5IX	Instant Fit CIC	36	10 - 80	4	4,33,990
Kit Silk C&G 3IX	Instant Fit CIC	32	10 - 80	2	3,14,990
Kit Silk C&G 2IX	Instant Fit CIC	24	10 - 80	2	1,99,990
Kit Silk C&G 1IX	Instant Fit CIC	16	10 - 80	2	1,39,990
KIT INSIO C&G CIC IX ^{New}					
Kit Insio C&G CIC 7IX ^{New}	CIC	48	10 - 110	4	7,69,990
Kit Insio C&G CIC 5IX ^{New}	CIC	36	10 - 110	4	4,79,990
Kit Insio C&G CIC 3IX ^{New}	CIC	32	10 - 110	2	3,46,990
Kit Insio C&G CIC 2IX ^{New}	CIC	24	10 - 110	2	2,04,990
Kit Insio C&G CIC 1IX ^{New}	CIC	16	10 - 110	2	1,44,990
INSIO C&G CIC IX ^{New}					
Insio C&G CIC 7IX ^{New}	CIC	48	10 - 110	4	3,89,990
Insio C&G CIC 5IX ^{New}	CIC	36	10 - 110	4	2,49,990
Insio C&G CIC 3IX ^{New}	CIC	32	10 - 110	2	1,79,990
Insio C&G CIC 2IX ^{New}	CIC	24	10 - 110	2	1,07,990
Insio C&G CIC 1IX ^{New}	CIC	16	10 - 110	2	77,990
INSIO IX					
Insio 7IX	ITC/ITE/CIC/IIC	48	10 - 100	4	3,25,990
Insio 5IX	ITC/ITE/CIC/IIC	36	10 - 100	4	1,94,990
Insio 3IX	ITC/ITE/CIC/IIC	32	10 - 100	2	1,15,990
Insio 2IX	ITC/ITE/CIC/IIC	24	10 - 100	2	74,990
Insio 1IX	ITC/ITE/CIC/IIC	16	10 - 100	2	51,990

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.

Product Name	Form Factor	Battery Type	Warranty (Years)	MRP (₹)
CROS / BiCROS				
Kit Active Pro CROS IX	Instant Fit ITE	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	2,12,990
Kit Active CROS IX	Instant Fit ITE	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,29,990
Kit Styletto CROS 7IX	Slim-RIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	4,12,990
Kit Styletto CROS 5IX	Slim-RIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	2,99,990
Kit Styletto CROS 3IX	Slim-RIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	2,24,990
Kit Styletto CROS 2IX	Slim-RIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,54,990
Kit Styletto CROS 1IX	Slim-RIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,24,990
Kit Silk C&G CROS 7IX	Instant Fit CIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	4,12,990
Kit Silk C&G CROS 5IX	Instant Fit CIC	Li-ion	4 yrs for Hearing Aid 1 yr for CROS	2,66,990
Kit Silk C&G CROS 3IX	Instant Fit CIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	2,07,990
Kit Silk C&G CROS 2IX	Instant Fit CIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,49,990
Kit Silk C&G CROS 1IX	Instant Fit CIC	Li-ion	2 yrs for Hearing Aid 1 yr for CROS	1,19,990
CROS / BiCROS (Pure C&G IX Transmitter)	RIC	Li-ion	1	35,990
CROS / BiCROS (Pure C&G BCT IX Transmitter)	RIC	Li-ion	1	35,990



Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.

Product Name		Warranty (Year)	MRP (₹)
CHARGERS			
Active IX Charger		1	11,990
Styletto IX Charger		1	11,990
Dry&Clean Charger Styletto		1	11,990
Dry&Clean Charger Pure		1	11,990
Silk Charger		1	11,990
Pure Portable Charger		1	7,990
Desktop Charger Pure		1	5,990
Multi-Charger		1	4,490
Desktop Charger Motion M/P/SP		1	5,990
Dry&Clean Charger Motion P/SP		1	11,990

Kit Includes 2 units of hearing aids and a standard charger.
 Note: The fitting range is completely dependent on the receiver matrix of the hearing aid



AX

Platform Features

Features	Information	Performance levels				
		7AX	5AX	3AX	2AX	1AX
Dynamic Soundscape Processing 2.0	Dynamic Soundscape Processing takes sound processing to a more precise level by including motion data for a more realistic sound adapted to everyday life.				...	...
Augmented Focus™	Separates the focus sound (speech) from surroundings, creating a clear contrast between the two, which makes it easier to follow conversations.	✓	✓	✓	✓	✓
Acoustic Sensors	The redesigned acoustic sensors can detect more variables of sound than ever before, providing more accurate information about the environment.	✓	✓	✓	✓	✓
e2e wireless 4.0	Control coupling, synchronization, binaural directionality, eWindscreen binaural, Spatial Speech Focus, CROS/ BiCROS.	✓	✓	✓	✓	✓
Motion Sensor	An integrated motion sensor detects whether the wearer is walking, running or sitting still. This motion data complements the acoustic data to optimize Dynamic Soundscape Processing 2.0.	✓	✓	✓		
OVP™ (Own Voice Processing)	Own Voice Processing (OVP™) utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression. Available for RIC/BTE devices with e2e wireless in binaural fittings.	✓	✓	✓		
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearing situations.			...	...	..
Signal processing / Gain & MPO	Channels/handles in Connexx. For Signia AX, channels are per processing chain (Focus/Surroundings).	48/20	32/16	24/12	16/8	16/8
Hearing Programs	Number of programs configurable in Connexx.	6	6	6	4	4
Extended dynamic range	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓
Speech and noise management	Reduces listening effort using sophisticated sound processing	✓	✓	✓	✓	✓
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	✓	✓	✓	✓	
Feedback cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓
HD Music	Preset programs for enhanced enjoyment of non-streamed music.	3	3	1	1	
eWindScreen	Reduces the annoyance of wind noise when outdoors.	✓	✓	✓	✓	
Extended bandwidth	12 kHz bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓				
EchoShield	Dedicated program for reverberant environments.	✓	✓			
Speech Quality	Speech Quality delivers improved speech understanding and communication ease in noise.			...	...	...
Binaural Directionality	Directionality provides scientifically proven enhanced speech understanding in noise. Available for bilateral fittings with directional microphones and e2e wireless link.	✓	✓	✓		

1)

In Universal program: Automatic mic beam to the right / left / rear directions.

2)

In Universal program: Automatic mic beam to the rear, surround mode for side speech.
In manual programs, microphone beam to the right / left / rear directions.

Features	Information	Performance levels				
		7AX	5AX	3AX	2AX	1AX
Wireless CROS/BiCROS	Requires CROS transmitters.	✓	✓	✓	✓	✓
Frequency compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓
Spatial SpeechFocus	Depending on the direction of the dominant speech source, Spatial SpeechFocus steers the directional beam to the left, right, or behind the wearer automatically in a car situation, or in a dedicated "Strolling" program. Available for bilateral fittings with directional microphones and e2e wireless link.	✓ ¹⁾	✓ ²⁾			
Wearer Interaction	Technically inclined wearers may appreciate having extra control over their hearing aids. Signia Augmented Xperience hearing aids offer these wearers a comprehensive range of adjustment options via the Signia App.					
Signia Assistant ¹	Signia Assistant offers wearers instant 24/7 support and troubleshooting around the clock. It provides detailed information in Connexx to address each client's fine-tuning needs more precisely	✓	✓	✓	✓	✓
Signia app (iOS and Android)	The new Signia app combines all functionalities of previous apps (touchControl, myControl, myHearing app). The functionality offered in the app is automatically customized depending on the connected hearing aid.	✓	✓	✓	✓	✓
Adaptive Streaming Volume	Speech intelligibility weighted adaptation of the streaming volume depending on ambient noise.	✓	✓	✓	✓	✓
Spatial Configurator	Spatial Configurator is part of the Signia app and allows the wearer to adjust the direction of the microphone beam.	✓	✓			
Direct streaming for iOS	Signia Augmented Xperience hearing aids can be connected directly to Apple devices for phone calls and audio streaming – without additional accessories.	✓	✓	✓	✓	✓
Direct streaming for Android	Direct audio streaming from ASHA compatible Android devices to Bluetooth enabled Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
StreamLine TV accessory	Stream high quality stereo sound directly from a TV to Bluetooth equipped Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
StreamLine Mic accessory	Stream high quality stereo sound directly from any Bluetooth device (Android devices, older Apple devices, other smartphones, PCs, laptops,etc.) to Bluetooth equipped Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
Tinnitus	A choice of tinnitus treatments based on amplified notch therapy or traditional noise therapy signals. It is even possible to combine the two approaches.	✓	✓	✓	✓	
Notched AmplificationTherapy	Innovative tinnitus therapy delivered via everyday amplification.	✓	✓	✓	✓	
Smart Optimizer and Data Logging	Smart Optimizer (found under Data Logging) provides tailored fine-tuning recommendations based on an analysis of the wearer's usage of thehearing aid.	✓	✓	✓	✓	✓
Acclimatization manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓
InSituGram	Allows quick and convenient screening audiograms and checking of thresholds while the hearing aids are worn.	✓	✓	✓	✓	✓
AutoFit	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets. The new IMC2 standard expands the range of REM equipment which supports AutoFit with Connexx and Noah 4.5.1 or higher.	✓	✓	✓	✓	✓
TeleCare	TeleCare provides remote services including additional tools to follow-up with patients. In addition, user engagement and autonomy can be increased via the Signia app.	✓	✓	✓	✓	✓
Remote Services	Includes TeleCare Dashboard, Remote Fine Tuning, CareChat (incl. text messages, phone calls and video calls), virtual follow-up.		✓	✓	✓	✓



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
AUGMENTED XPERIENCE					


KIT STYLETTO AX					
Kit Styletetto 7AX	Slim-RIC	48	10 - 110	4	6,99,990
Kit Styletetto 5AX	Slim-RIC	32	10 - 110	4	4,68,990
Kit Styletetto 3AX	Slim-RIC	24	10 - 110	2	3,42,990
Kit Styletetto 2AX	Slim-RIC	16	10 - 110	2	1,99,990
Kit Styletetto 1AX	Slim-RIC	16	10 - 110	2	1,39,990

KIT PURE C&G AX / KIT PURE C&G T AX					
Kit Pure C&G 7AX / T 7AX	RIC	48	10 - 115	4	7,45,990
Kit Pure C&G 5AX / T 5AX	RIC	32	10 - 115	4	4,64,990
Kit Pure C&G 3AX / T 3AX	RIC	24	10 - 115	2	3,35,990
Kit Pure C&G 2AX / T 2AX	RIC	16	10 - 115	2	1,89,990
Kit Pure C&G 1AX / T 1AX	RIC	16	10 - 115	2	1,34,990

PURE C&G AX / PURE C&G T AX					
Pure C&G 7AX / T 7AX	RIC	48	10 - 115	4	3,80,990
Pure C&G 5AX / T 5AX	RIC	32	10 - 115	4	2,39,990
Pure C&G 3AX / T 3AX	RIC	24	10 - 115	2	1,75,990
Pure C&G 2AX / T 2AX	RIC	16	10 - 115	2	99,990
Pure C&G 1AX / T 1AX	RIC	16	10 - 115	2	72,990

PURE 312 AX					
Pure 312 7AX	RIC	48	10 - 115	4	3,30,990
Pure 312 5AX	RIC	32	10 - 115	4	2,09,990
Pure 312 3AX	RIC	24	10 - 115	2	1,49,990
Pure 312 2AX	RIC	16	10 - 115	2	74,990
Pure 312 1AX	RIC	16	10 - 115	2	52,990

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
AUGMENTED XPERIENCE					

KIT INSIO C&G AX					
Kit Insio C&G 7AX	ITC / ITE	48	10 - 110	4	7,45,990
Kit Insio C&G 5AX	ITC / ITE	32	10 - 110	4	4,64,990
Kit Insio C&G 3AX	ITC / ITE	24	10 - 110	2	3,35,990
Kit Insio C&G 2AX	ITC / ITE	16	10 - 110	2	1,89,990
Kit Insio C&G 1AX	ITC / ITE	16	10 - 110	2	1,34,990

INSIO C&G AX					
Insio C&G 7AX	ITC / ITE	48	10 - 110	4	3,80,990
Insio C&G 5AX	ITC / ITE	32	10 - 110	4	2,39,990
Insio C&G 3AX	ITC / ITE	24	10 - 110	2	1,75,990
Insio C&G 2AX	ITC / ITE	16	10 - 110	2	99,990
Insio C&G 1AX	ITC / ITE	16	10 - 110	2	72,990

INTUIS 4					
Intuis 4.7 (SP/P/M)	BTE	48	10 - 120	4	2,15,990
Intuis 4.5 (SP/P/M)	BTE	32	10 - 120	4	1,36,990
Intuis 4.3 (SP/P/M)	BTE	24	10 - 120	2	79,990
Intuis 4.2 (SP/P/M)	BTE	16	10 - 120	2	43,990
Intuis 4.1 (SP/P/M)	BTE	16	10 - 120	2	37,990
Intuis 4.0 (SP/P/M)**	BTE	12	10 - 120	2	31,990

Note:
The fitting range is completely dependent on the receiver matrix of the hearing aid.
** Direct streaming option not available in Intuis 4.0.



Product Name	Form Factor	Battery Type	Warranty (Years)	MRP (₹)
CROS / BiCROS				
Kit Styletto CROS 7AX	Slim-RIC	Li-ion	4 years for Hearing Aid 1 year for CROS	3,99,990
Kit Styletto CROS 5AX	Slim-RIC	Li-ion	4 years for Hearing Aid 1 year for CROS	2,84,990
Kit Styletto CROS 3AX	Slim-RIC	Li-ion	2 years for Hearing Aid 1 year for CROS	2,21,990
Kit Styletto CROS 2AX	Slim-RIC	Li-ion	2 years for Hearing Aid 1 year for CROS	1,49,990
Kit Styletto CROS 1AX	Slim-RIC	Li-ion	2 years for Hearing Aid 1 year for CROS	1,19,990

Product Name		Warranty (Year)	MRP (₹)
CHARGERS			
Styletto AX Charger		1	11,990
Dry&Clean Charger Styletto		1	11,990
Desktop Charger Insio AX		1	11,990
Multi-Charger		1	4,490



Orion Charge&Go feature overview

Features	Performance levels			
	200	100	75	50
Frequency channels / compression controls	16/8	16/8	12/6	8/4
Number of programs	4	4	4	4
48-channel first fit	●	●	●	●
Adaptive compression	●	●	●	●
Extended dynamic range	●	●	●	●
Signia app	●	●	●	●
miniPocket compatability	●	●	●	●
InSituGram	●	●	●	●
AutoFit	●	●	●	●
Smart optimizer and Data logging	●	●	●	●
Feedback cancellation	●	●	●	●
Speech and noise management	●	●	●	●
Microphone directionality	●	●	●	●
Automatic directionality	●	●	●	-
Adaptive directionality	●	-	-	-
TruEar	●	●	-	-
SoundBalance	●	●	-	-
Frequency compression	●	●	-	-
Automatic acclimitization manager	●	●	-	-
e2e wireless	●	●	-	-
User control coupling via e2e	●	●	-	-
SoundSmoothing™	●	-	-	-
eWindScreen™	●	-	-	-
Tinnitus static therapy signal	●	-	-	-

Effortless to use

- Rechargeability for all
Your clients can say goodbye to the hassle of small batteries and hello to rechargeability.
- Wireless fitting
Cumbersome cables are a thing of the past with wireless technology that streamlines the fitting process.
- First fit precision
Our 48 channel first fit calculation . ensures a precise and personalized fit from the start, reducing the need for fine-tuning and optimizing the client experience right from the start.

Elevated design

- Contemporary design
With an attractive premium finish, Orion Charge&Go's design is as sophisticated as its technology.
- IP68 rating
Certified to resist dust and sweat, Orion Charge&Go provides dependable performance in all conditions.

Enhanced sound processing

- Clarity in quiet
Our TruEar technology emulates the effect of the natural ear shape by prioritizing speech frequencies from the front for optimal hearing support in quiet environments.
- Focus in noise
Adaptive directionality zeros in on the speech your client wants to hear in noisy environments.
- Versatility across soundscapes
Adaptive compression and an extended input dynamic range work in tandem to adjust amplification to the natural flow of sounds, ensuring impeccable sound quality everywhere.



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
ORION C&G P/SP BTE					
Orion C&G P/SP 200	BTE	16	10 - 120	2	42,990
Orion C&G P/SP 100	BTE	16	10 - 120	2	40,990
Orion C&G P/SP 75	BTE	12	10 - 120	2	35,990
Orion C&G P/SP 50	BTE	8	10 - 120	2	32,990
ORION C&G P/SP T BTE					
Orion C&G P/SP T 200	BTE	16	10 - 120	2	42,990
Orion C&G P/SP T 100	BTE	16	10 - 120	2	40,990
ORION C&G RIC					
Orion C&G 200	RIC	16	10 - 115	2	42,990
Orion C&G 100	RIC	16	10 - 115	2	40,990
Orion C&G 75	RIC	12	10 - 115	2	35,990
Orion C&G 50	RIC	8	10 - 115	2	32,990

Product Name		Warranty (Year)	MRP (₹)
CHARGER			
Multi-Charger		1	4,490



Note:
Charger for Orion C&G will have to be purchased separately.
The fitting range is completely dependent on the receiver matrix of the hearing aid.

Connect with your expectations

With Sirion Connect, we're once again redefining what's possible in the essentials category, bringing seamless connectivity, exceptional design and advanced technology to more people than ever before.

Effortless to use

- **Connectivity for more**
So more of your clients can stream music.
- **Tailored control**
Your clients can adjust the settings to their needs with the user-friendly Signia app.
- **Wireless fitting**
For a streamlined fitting process without the hassle of cables.
- **First fit precision**
Our 48 channel first fit approach ensures a precise and personalized fit from the start, reducing the need for fine-tuning.

Elevated design

- **Contemporary design**
Sirion Connect combines a premium finish with easy-to-use push buttons.
- **IP68 rating**
Featuring an advanced nanocoating that repels water, Sirion Connect is IP68-certified to resist dust and sweat, even in the most challenging environments.

Enhanced sound processing

- **Split processing**
We are bringing our advanced sound processing technology to the essentials category. By handling speech and background sound separately, split processing makes speech crisp while keeping ambient sounds comfortable. Proven to enhance speech understanding in busy environments, it lets wearers stay engaged with ease.
- **Clarity in quiet**
Split processing technology prioritizes speech frequencies from the front for crystal clear conversations.
- **Focus in noise**
As it gets louder, split processing ensures that background sound is more aggressively reduced while separately ensuring speech is amplified for greater contrast.
- **Versatility across soundscapes**
Split processing stays on, constantly adapting to the environment to ensure speech can always be heard separately from background sounds.

Connect with your expectations



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
SIRION CONNECT					
Sirion Connect 75	RIC	12	10 - 115	2	32,990
Sirion Connect P 75	BTE	12	10 - 115	2	32,990
Sirion Connect SP 75	BTE	12	10 - 120	2	32,990

Connect with more

Connect with your expectations



Sirion Connect

Effortless. Elevated. Enhanced.



Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
ESSENTIAL					
Intuis 3 Custom	ITE/ITC/CIC	12	10 - 110	2	31,990
Intuis 3 312	RIC	12	10 - 115	2	26,990
Intuis 3 SP	BTE	12	20 - 120	2	26,990
Intuis 3 P	BTE	12	10 - 110	2	26,990
Intuis 3 M	BTE	12	10 - 90	2	26,990
Intuis 3 S	BTE	12	0 - 80	2	26,990
Intuis 3 Click	ITC/CIC	12	10 - 90	2	26,990
BASIC					
Prompt Custom	ITC/CIC	8	10 - 110	2	26,990
Prompt SP	BTE	8	20 - 120	2	20,990
Prompt P	BTE	8	10 - 110	2	20,990
Prompt S	BTE	8	0 - 80	2	20,990
Prompt Click	ITC/CIC	8	10 - 90	2	20,990
Run Click	ITC/CIC	8	10 - 80	2	18,990
Run SP	BTE	8	60 - 120	2	18,990
Run P	BTE	8	20 - 100	2	18,990
Fun SP	BTE	6	60 - 120	2	13,990
Fun P	BTE	6	20 - 100	2	13,990
Fast P	BTE	4	20 - 100	2	10,990

Product Name	Form Factor	Channels	Fitting Range (dB)	Warranty (Years)	MRP (₹)
OTHER HEARING AIDS					
Kit Motion C&G 7X / P 7X / SP 7X	BTE	48/48/34	10 - 120	4	6,50,990
Kit Motion C&G 5X / P 5X / SP 5X	BTE	32/32/32	10 - 120	4	4,09,990
Kit Motion C&G 3X / P 3X / SP 3X	BTE	24/24/24	10 - 120	2	2,30,990
Kit Motion C&G 2X / P 2X / SP 2X	BTE	16/16/16	10 - 120	2	1,36,990
Kit Motion C&G 1X / P 1X / SP 1X	BTE	16/16/16	10 - 120	2	1,04,990
Motion C&G 7X / P 7X / SP 7X	BTE	48/48/34	10 - 120	4	3,32,990
Motion C&G 5X / P 5X / SP 5X	BTE	32/32/32	10 - 120	4	2,12,990
Motion C&G 3X / P 3X / SP 3X	BTE	24/24/24	10 - 120	2	1,18,990
Motion C&G 2X / P 2X / SP 2X	BTE	16/16/16	10 - 120	2	75,990
Motion C&G 1X / P 1X / SP 1X	BTE	16/16/16	10 - 120	2	59,990
Silk 7X	Instant Fit CIC	48	10 - 80	4	3,14,990
Silk 5X	Instant Fit CIC	32	10 - 80	4	1,88,990
Silk 3X	Instant Fit CIC	24	10 - 80	2	1,36,990
Silk 2X	Instant Fit CIC	16	10 - 80	2	62,990
Silk 1X	Instant Fit CIC	16	10 - 80	2	46,990
Insio 7NX	ITE/ITC/CIC/IIC	48	10 - 110	4	3,14,990
Insio 5NX	ITE/ITC/CIC/IIC	32	10 - 110	4	1,89,990
Insio 3NX	ITE/ITC/CIC/IIC	24	10 - 110	2	1,09,990
Insio 2NX	ITE/ITC/CIC/IIC	16	10 - 110	2	59,990
Insio 1NX	ITE/ITC/CIC/IIC	16	10 - 110	2	44,990
Product Name	Form Factor	Battery Type	Warranty (Year)		MRP (₹)
CROS / BiCROS					
CROS / BiCROS (Silk X Transmitter)	Instant Fit CIC	10	1		29,490
CROS / BiCROS (Pure C&G AX Transmitter)	RIC	Li-on	1		35,990
CROS / BiCROS (Pure 312 AX Transmitter)	RIC	312	1		29,490

Accessories for Hearing Aids			
Product	MRP (₹)	Product	MRP (₹)
Remote Control and Accessories			
TV Sound	17,700	EarWear 2.0	
StreamLine Mic	22,000	Domes	
StreamLine TV	19,470	Click Domes Open (8mm, 10mm) (6pcs/blister)	932
miniPocket	7,670	Click Domes Closed (8mm, 10mm) (6pcs/blister)	932
Hearing Aid Zinc Air Batteries		Click Domes Double (8mm/10mm,10/12 mm) (6pcs/blister)	932
Signia # 675 Battery (Mercury Free)	51	Click Domes Open (8mm/10mm) (6pcs/blister)	932
Signia # 13 Battery (Mercury Free)	51	Sleeves	
Signia # 312 Battery (Mercury Free)	51	Sleeves Vented (XS,S,M,L) (6pcs/blister)	1,168
Signia # 10 Battery (Mercury Free)	51	Sleeves Closed (XS,S,M,L) (6pcs/blister)	1,168
In-The-Ear (ITE) Fitting Accessories		miniReceiver 2.0	
Quick Guard Kit	530	S miniReceiver 2.0 (Sizes 0 - 3) Left / Right	5,890
Receiver-In-Canal (RIC) Fitting Accessories		M miniReceiver 2.0 (Sizes 0 - 4) Left / Right	5,890
Receiver Mold	1,168	P miniReceiver 2.0 (Sizes 0 - 3) Left / Right	5,890
miniReceiver Waxguard Set (10 pcs)	650	HP miniReceiver 2.0 (Sizes 0 - 3) Left / Right	5,890
Behind-The-Ear (BTE) Fitting Accessories		EarWear 3.0	
Comfort Coupling Set (for Prompt, Run, Fun)	236	Sleeves 3.0	
Simple Coupling Set (for Fast)	177	Sleeves 3.0 Vented (XS, S, M, L) (6pcs/blister)	1,190
BTE Eartip Assembly (Eartip + Elbow + Tube)	118	Sleeves 3.0 Closed (XS, S, M, L) (6pcs/blister)	1,190
BTE Eartip Joint Elbow (L-Bend)	41	Sleeves 3.0 Power (XS, S, M, L) (6pcs/blister)	1,190
BTE Eartip (Any Size)	59	Eartips 3.0	
BTE Hucki Holder (One Pair)	148	Eartips 3.0 Open (5mm, 7mm, 10mm) (6pcs/blister)	790
Colour Exchange Kit	1,522	Eartips 3.0 Closed (8mm, 12mm) (6pcs/blister)	790
Miscellaneous Items		ThinTube 3.0	
UV Hard Pocket / BTE Earmold	900	Set Thintube 3.0 (S0, S1, S2, S3, S4) (Set of Left & Right)	490
Silicon Soft Pocket / BTE Earmold	1,100	Set Thintube 3.0 P (S0, S1, S2, S3, S4) (Set of Left & Right)	490
Custom Swim Plug	1,100	miniReceiver 3.0	
Drying Kit	150	S miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890
Cleaning Kit	590	M miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890
CROS / BiCROS Transmitter		P miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890
Styletto IX Transmitter	39,990	HP miniReceiver 3.0 (Sizes 0 - 4) Left / Right	5,890
Silk C&G IX Transmitter	39,990		
Active IX Transmitter	39,990		
Styletto AX Transmitter	39,990		

App, accessory, and charger compatibility

This overview shows you which apps and accessories work with our hearing aids. Depending on the solution your customers choose, they can even enhance their hearing experience and convenience with this well-matching portfolio.

Hearing aid	Signia app	CROS/ BiCROS	Streamers and transmitters		Remote Controls	Portable Charger	Inductive Charger II	Desktop Charger	Dry&Clean Charger	Multi Charger
										
Styletto IX	●	● ^A	●	●	●	●	-	-	●	-
Pure Charge&Go IX family	●	● ^B	●	● ^F	●	●	-	●	●	●
Active Pro IX	●	● ^C	●	●	●	●	-	-	-	-
Motion Charge&Go IX	●	● ^D	●	●	●	-	-	●	●	●
Silk Charge&Go IX	●	● ^E	-	-	●	●	-	-	-	-
Insio IX family	●	● ^E	-	-	●	● ^G	-	-	-	-
Styletto AX	●	● ^H	●	●	●	●	-	-	●	-
Insio Charge&Go AX	●	● ^{I,J}	●	●	●	-	-	●	-	-
Pure Charge&Go AX / T AX	●	● ^{I,J}	●	●	●	●	-	●	●	●
Pure 312 AX	●	● ^{I,J}	●	●	●	-	-	-	-	-
Motion Charge&Go X	●	● ^K	●	●	●	-	●	-	-	-
Motion Charge&Go P X	●	● ^K	●	●	●	-	-	●	●	●
Motion Charge&Go SP X	●	● ^K	●	●	●	-	-	●	●	●
Silk X	●	● ^L	-	-	●	-	-	-	-	-
Pure 312 X	●	● ^K	●	●	●	-	-	-	-	-
Sirion Connect	●	-	●	●	●	-	-	-	-	●
Orion Charge&Go	●	-	-	-	●	-	-	-	-	●
Intuis 4	●	● ^{I,J,M}	● ^M	● ^M	●	-	-	-	-	-
Intuis 3 RIC312	●	-	-	-	●	-	-	-	-	-
Intuis 3	●	-	-	-	●	-	-	-	-	-
Prompt	●	-	-	-	●	-	-	-	-	-
Run	●	-	-	-	-	-	-	-	-	-

- A) Compatible with CROS Styletto IX
- B) Pure Charge&Go IX / T IX are compatible with CROS Pure Charge&Go IX; Pure Charge&Go BCT IX is compatible with CROS Pure Charge&Go BCT IX
- C) Compatible with CROS Active IX
- D) Compatible with CROS Pure Charge&Go IX

- E) Compatible with CROS Silk Charge&Go IX
- F) Not compatible with Pure Charge&Go BCT IX
- G) Only for Insio Charge&Go CIC IX
- H) Compatible with CROS Styletto AX
- I) Compatible with CROS Pure Charge&Go AX
- J) Compatible with CROS Pure 312 AX
- K) Compatible with CROS Pure 312 X
- L) Compatible with CROS Silk X

- M) Excluding PL 4.0
- N) Compatible with Bluetooth-enabled Insio Nx ITC/ITE